

# HM-NEB-ULTRASONIC



**HM-NEB-ULTRASONIC** The ultrasonic nebulizer can nebulize the liquid medication into tiny mist particles of uniform size, which can be inhaled by patients to treat the diseases. This has been proven clinically effective for treatment of coryza, pharyngitis, bronchitis, etc.

## Ultrasonic Nebulizer

**Features:** new generation of upgraded product. Durable. Stretchable Double tubes can produce a dense and tiny-particle mist. (Major components such as circuit board etc are connected as socket connectors, more convenient for repair and replacement, and equipped with a package for outside use)

### **Principal technical data:**

Power Voltage: AC220V  $\pm$  10%

Ultrasonic frequency: 1.7MHz  $\pm$  10%

Power:  $\leq$ 60W

Maximum nebulizing rate: 4ml/min

Timing range: 0-60 min no level for adjustment

Continuous loading time:  $\geq$  4h

Max medication cup capacity: 150mL

### **HILLMED CORPORATION**

7215 N.W. 46th St. Miami, FL 33166  
Ph: (305) 594-7474 Fax: (305) 477-0699

[www.hillusa.com](http://www.hillusa.com)

